




A NECESSITY FOR HOME. A SAFE AND CONVENIENT ELECTROMOTIVE EAR-CLEANER, EAR CLEAN.



## The first and newest invention in the world, **Ear Cleaner**

### Why do people love our ear cleaner?

The ear clean is a safe & convenient ear cleaner which helps maintain & prevents ear from all types of ear diseases. By sucking ear-waxes and humidity out of the ear, its unique design by using a revolving inhalation screw to protect against any inflammation that may cause increasing virus to the ears.

- 1 When water or foreign material like sand get into your ears, after showers or swimming. The sucking action will pull them out.
- 2 When small insects such as flies, ants, mosquito, etc.. Irritate the ear, no problem with the ear clean.
- 3 The ear clean can make a great gift for friends & family.



#### ► At Home

You may hesitate to clean inside your children's ears because their small ears have soft skin. The nontoxic silicon tube will cleanse inside your children's ears without any pain. An adult can cleanse his own ear without assistance.



#### ► At Barber Shop or Beauty Salon

These days, beauty salon offer skin treatment services. A clean face with dirty ears may not satisfy your customer. The ear clean can solve the problem by cleansing the ear. You can offer better service to your customer by using ear clean.



#### ► At Hospital and Office

Until now, doctors have gathered up the dissolved Ear-waxes of patients by using chemical which may irritate the ear. ear clean is a safe and painless way to provide a healthy ear.



# The Ear Clean, which was designed to match the structure of human body, will give you great satisfaction.

The Ear Clean, which we invented by putting our years of engineering experience into detailed parts, is the first and newest device in the world, which is gaining public favor after being imported to advance countries like Japan.

## Nontoxic Silicon Tube

The nontoxic silicon tube was designed engineered for the ear, which will cleanse your children's ears without any pain.

## Inhalation Screw

The inhalation screw which is connected to the motor was designed to suck up separated ear-wax by the silicon tube without hurting your eardrum.

## Ear-Wax Box

This gathers the separated ear-waxes so you can clean this up easily with attached brush after opening its cap.

## Filter

This prevents small ear-waxes from entering into the motor expanding the life span of the Ear Clean. It can be separated from the main body easily to be cleansed with water.

## ■ Specification

Voltage	DC 3V	LR6 AA Size 1.5V × 2EA
Electric Current	0.5 A	
Sound-level	50dB under	
R.P.M	14,600	Fan Load
Function	ON / OFF	2 Level switch
Weights	107g	NET
Materials	ABS+Silicone	
Product Size	100(D) × 59(W) × 130(H)	1C/T : 40ea
Cin f Size	109(D) × 65(W) × 139(H)	Box type
	155(D) × 62(W) × 250(H)	Buster type

**Recommended Battery**  
Normal AA 1.5V battery × 2

